vivocare  
health enjoyment

## Model: Fondly

**EN** Electronic Baby Scale

**TR** Elektronik Bebek Baskülü

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**IT** Bilancia pesaneonati elettronica

**DE** Elektronische Babywaage

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**EN** Instruction Manual

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## Features

- 1.4 High precision strain gauge sensor.
2. Auto on and auto off.
3. Tare function.
4. Auto lock/hold function.
5. Low battery and over load indication.
6. Unit conversion switch.

## Product Specification

Product Name: Electronic Baby Scale

Model: Fondly

LCD size: 60 x 20mm

Power supply: 2 xAAA battery

Measuring rang:100g-20kg

Precision:10g

## Operating Instructions

### 1. Battery and Battery Installation

This product adopts the environmental protection design of low power dissipation and has contained two pcs of AAA (No.7) batteries when leaving the plant.

Please tear off the thin insulation film of the batteries and install them directly,then you can use the products.Please change the batteries when the low voltage prompt shows.

- A. Open the cover of the battery compartment located at the bottom of the electronic scale.
- B. When installing new batteries,insert one end of the battery in the compartment first, then press the other end into position ensuring correct polarity.

## Instruction Diagram

## Button Functions

A. On/Off Button: Press this button to turn the scale On or Off.

B. Tare Button: This will allow you to add blankets or other items and then subtract their weight before adding your child.

C. The scale has various units from which to choose. Press the "UNIT" Button to choose.

### Weighing Description

1. Place the scale on a flat, hard surface.
2. Press "ON/OFF" button to turn the scale on, the display shows "Hi".
3. After LCD shows "0.00kg", it could be able to measure weight.
4. Press the "UNIT" Button to switch the displayed units between gram, kilograms, pound, stone and pound.
5. Put infant on the scale platform. The LCD will display the baby's weight.
6. Press "ON/OFF" key to turn off the scale or scale will turn off itself after 30 seconds.

### Auto lock function

When scale has the weight of the stuff, do not move the stuff off the platform, it will auto lock after Three flickers.

### Overload Indication

If the weight is over loaded "max capacity" LCD display will show " 0-Ld".

### Low Battery Indication

When LCD display shows "LO", Please replace the new battery with correct polarity.

### Important Safeguards

1. Place the scale on a flat surface (i.e a table) before use.
2. Do not subject the scale to very hot or cold environments.
3. Ensure that the scale does not fall as this could lead to damage as it consists of high precision electronics.
4. Use only a soft cleaning cloth to clean the product. Do not use any chemical solutions and keep the scale away from water as it is not waterproof.
5. This product is for household use only.
6. Please remove the batteries if the scale is not going to be used for a long time.
7. Please change the batteries if "LO" displays on the LCD.
8. If the LCD shows "0-Ld" please immediately remove the object from the scale as this shows that a overload has occurred.
9. Avoid shaking or hitting the scale please.
10. When in storage, please do not place anything on the scale.
11. Not for legal trade.

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