
promedic by comen

promedic by comen	C60
第任儿 → ¹⁹⁶⁹⁻¹²⁻³¹ 11 x1 0前 11 x1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	ECG PACE A 120113 512 0.00 A 50 971 0.00 A 51 0.00 0.20 9000 0.20 10 0 0
	NIBP 19:00 moldy 120 90 130 ^{MS} 198 ≇35
	SP02 PR 60 138 E 98 50 X RESP
●●●● ●●●	20 20 。 NBP湖堡 》 主来单 主屏幕

C60 Neonatal Monitor





PROMEDIC by COMEN C60 is CE certified, using the world advanced technology, innovate and develop the world leading C60 neonatal monitor. C60 equips with 8.4 high-definition color TFT screen, providing the touch screen and information hand-w-riting input function, embody the user-friendly design with special care to all users. According to the character of neonatal cardiovascular System, PROMEDIC by COMEN develops ExNeo® Neonate ECG technology, Adap-DSP® NIBP measurement technique, caring for the life with precision and profession attitude. Also, world famous OxiMax® neonatal SpO₂ technology, SatSecondsTM intellectualized alarming system, Specialized neonatal suffocation self-saving function equipped with oxygen density monitoring function; special newborn accessories are applied to the model to provide more care to the neonates.C60 model fills up the blank of the Chinese Neonatal Monitoring field, giving more energy to the world neonatal monitoring field.

C60 Neonatal Monitor







Various mounting solutions





- USB socket: extend memory capacity, convenient for computer data saving and file recording
- Net Port: Maximum 128units bedside monitors to connect the central monitoring system, wired or wireless network supportable

C60 Neonatal Monitor Protection




- C60 enhance the anti-interfere capability of the NIBP measurement to ensure the data accuracy in neonate movement status
- According to the low blood and low perfusion character of neonate, PROMEDIC by COMEN use Adap-DSP[®] system to provide more accurate measurement
- Single neonate measurement mode, soft & hard over-pressure protection



- OxyCRG: Display the interactive relation between heart rate, respiration and oxygen on the same screen, convenient for observing the clinic change of neonate
- ST analysis, Arr. analysis, Drug-dose calculation.etc

SatSeconds=Desaturation PointsXSeconds



CE qualified special neonatal wrapping SpO₂ sensor, anti-avulsion, avoid cross infection



- Lightweight and compact, 3.5KG
- Li-ion battery, 4 hours working time





- Apnea Self saving monitoring : Unique neonatal Apnea Self-saving monitoring, solve the suffocation symptom brought by Apnea neonatorum, gives superier care to the neonate.
- Oxy gen concentration monitoring function : Medical stuff can command and acquire change of incubator oxygen density through monitoring oxygen density with oxygen density sensor.



Share the human touch brought by versatile Neonatal Monitor, when

the life begins...

Standard Configuration:

3- lead ECG , HR , Masimo ${\rm SpO}_2$, PR , NIBP , RESP , Dual TEMP , bedrail hook, Vibro acoustic simulator , Oxygen concentration sensor

Optional Configuration: 5-lead ECG , Nellcor SpO₂ , Dual IBP , ETCO₂ , printer





İthalatçı/Distributor: Trimpeks İth.İhr.Tur.ve Tic.A.Ş. İstanbul-Türkiye, Tel +90 212 319 50 00 Fax +90 212 319 50 50, www.trimpeks.com